

FJR0021 - Olympus IX-81 Chassis Commands

Last updated 1 October 2009.

Introduction

The Olympus IX-81 is a fully automated microscope. If you are buying one, you should get it with the ZDC hardware autofocus system. It, or Nikon's Perfect Focus system, or an equivalent piece of equipment, should be standard on all microscopes with motorized objective position.

The core of the microscope is the chassis, which includes the entire optical path, plus the transmitted light source. Olympus refuses to release the specifications to communicate with the microscope chassis without a non-disclosure agreement. This is deeply unethical, but I have reverse engineered most of the command set to atone for their sin.

Things to keep in mind:

- When talking about the z-axis of a microscope, use "near" and "far" instead of "up" and "down." "Nearer" always means the objective ends closer to the sample; "farther" means the objective ends farther away. On an inverted microscope, "near" is up and "far" is down; on an upright microscope it is exactly the reverse. Better to use "near" and "far" to avoid confusion.
- You can always get the current state of the system by sending the command you would use to change that state followed by ?. For example, to get the current objective position, send 10B?. The microscope returns 10B 3, say, if the current objective is position 3 on the nosepiece.
- The microscope only understands positive integers, no negative numbers, no floating point. All distances are sent as positive integers measured in hundredths of a micron. All voltages are sent as tenths of a volt. Where negative numbers are needed, such as to specify relative motion, an extra argument is used to tell the microscope the sign of the number.

Communicating with the microscope

The IX-81 chassis is connected to a computer by an RS-232 serial cable. You need the right settings for the RS-232 communication, as reported by PortMon, are

```
Rate: 19200
StopBits: ERROR Parity: EVEN WordLength: 8
EOF:0 ERR:0 BRK:0 EVT:0 XON:0 XOFF:0
Shake:1 Replace:40 XonLimit:1 XoffLimit:1
RI:10 RM:1 RC:10 WM:100 WC:500
Mask: RXCHAR CTS
```

All commands are ASCII text terminated by CR LF (0x0d 0x0a, the DOS newline). The microscope responds in the same way. A command consists of a command name followed by a space then a series of comma separated arguments containing no spaces, or a command name followed directly by a question mark (?) to query the current state of what that command would change. For example, 1LMPSW ON turns on

the transmitted light lamp, `2MOV N,50,1,300000,49` issues a command to move the objective along the z-axis, and `1LMPSW?` asks for the current state of the transmitted light lamp.

When you query the state, the microscope returns the command, followed by a space, followed by the state. So `1LMPSW?` returns `1LMPSW ON` if the lamp is on or `1LMPSW OFF` if it is off.

When you issue a command to change the state, the microscope returns the command, followed by a space, followed by a plus sign (+). So `1LMPSW ON` returns `1LMPSW +` if it succeeds. If a command fails, the microscope responds with an x in place of the +.

You can continue to issue commands even while a current command is processing. You know a command is complete when you receive the response with that command's name and + or x.

The microscope ignores all lines of text it receives that don't begin with 1 or 2. Invalid commands that do start with 1 or 2 produce `1x` and `2x` as a response, respectively. For example,

```
1rubbish -> 1x
2rubbish -> 2x
```

Configuring the hardware

You can always use commands to get the state of the microscope, but if you want to change the state, you have to log in:

```
1LOG IN
2LOG IN
```

When you log in, it blocks the controls on the microscope unless you explicitly reenable them. This is why there are two separate commands. `1LOG` enables all commands beginning with 1, which control selecting condensers and objectives, turning lamps on and off, and opening and closing shutters. `2LOG` controls z-motion and autofocus. Having them separate means you can use the microscope's focus wheel as if it were a manual scope, while still controlling the choice of objective from the computer.

To reenable the focus wheel on the microscope after `2LOG IN`, use

```
2JOG ON
2JOGSNS 10
2joglmt ON
```

`2JOG` enables or disables the wheel. `2JOGSNS` sets its sensitivity, and `2joglmt` stops the wheel motion from going beyond the near and far limits of the objective position (see `2FARLMT` and `2NEARLMT`).

The `1LOG IN` line lets you affect the state of all commands beginning with 1; the `2LOG IN` line does the same for command beginning with 2. This is actually fairly reasonable. All commands to control lamps, condensers, objectives, and shutters begin with 1; z-axis motion of the objective and autofocus begin with 2.

Using the buttons on the chassis

The Olympus IX-81 has a number of buttons for focusing up and down, controlling the lamp, switching


between coarse and fine focus, and often a keypad with a number of buttons. When you have sent the login commands, none of these function as they did. You have to handle their events yourself.

First, send

1SW ON

which will make the IX-81 send a string everytime one of the buttons is pressed or released. It uses the same string for any button release, so you cannot work with any combinations of buttons, only one at a time. When any button is released, the IX-81 sends 1SW 0.

When a button is pressed it sends 1SW *n* where *n* is an integer indicating the button. The buttons are (written in scientific notation -- for actual use, write them out with no commas, so 1e4 will appear as 10000)





If you do not have precisely the keypad shown above, this is no great problem. Simply run HyperTerminal or Terminate, log in and turn on 1SW ON, then start hitting the buttons and write down the numbers which are sent to the terminal.

To disable the button commands, send 1SW OFF.

Objectives, condensers, etc.

This section is organized along the lightpath, beginning at the lamp. Only light sources such as the transmitted light lamp which are mounted directly on the chassis are controlled by this command set. Most mercury or xenon lamps for fluorescence microscopy are external. If you do have a second internal light source, the command to change between light sources is 1LMPSEL. To select the transmitted light source, issue

```
1LMPSEL DIA
```

I have no other light sources on my chassis, and so no idea what the relevant commands might be.

You can turn lamps on and off with

```
1LMPSW ON
1LMPSW OFF
```

Its intensity is controlled by the supplied voltage. For the transmitted light lamp the voltage ranges from 0V to 12V. The chassis accepts voltages given in tenths of volts, so setting the intensity to 0V, to 5.6V, and to 12V, would be

```
1LMP 0
1LMP 56
1LMP 12
```

The chassis's control panel has two shutter ports, labelled SHA1 and SHA2. Generally SHA1 is connected to the transmitted light shutter, and SHA2 to the fluorescence shutter. If you have additional light sources and shutters, they will be controlled directly, not through the chassis. To open and close shutter 1, then

shutter 2, issue

```
1SHUT1 OUT
1SHUT1 IN
1SHUT2 OUT
1SHUT2 IN
```

Mnemonically, a shutter is closed if it has put something in the light path, and open if it has taken it out. There is a third command to control the external shutter for a fluorescent light source: 1LED. It accepts one argument, either 1 or 2. Which one opens or closes the shutter depends on whether your particular external shutter is open or closed when it receives a low voltage on its control wire (2 is low voltage; 1 is high voltage).

If you have a motorized condensor, 1CD # selected item # in the condensor wheel.

The positions on your objective nosepiece are numbered 1, 2, etc. To switch to the objective in position 3, issue

```
1OB 3
```

Beyond the objective is your filter cube. Again, the positions in the wheel containing these cubes are labelled 1, 2, etc. To select the filter in position 2, issue

```
1MU 2
```

Beyond the filter cube is the prism which directs the light either to the camera port or to the eyepiece. 1PRISM 1 shifts to the eyepiece, 1PRISM 2 to the camera. If you have a lower back port for an additional camera, 1BPORT controls this, but I don't have one and can't test it. For safety reasons, you should always close the external fluorescence shutter with 1LED when you use 1PRISM to switch to the eyepiece.

Motion

The only motion controlled by the chassis is the z-position of the objective. Before you start moving it, you should set the near and far limits, such as with

```
2FARLMT 10
2NEARLMT 3000000
```

To get the current position, use 2POS?. The microscope returns something like 2POS 539031, where 539,031 is the current z-position in hundredths of a micron from its farthest position.

To actually move the objective, use 2MOV.

```
2MOV F,2500,1,300000,49
2MOV N,300,1,300000,49
2MOV d,539031,1,300000,49
```

The first argument is one of N, F, and d, and it changes the meaning of the second argument. N and F make the second argument a distance to move relative to the current position, N to move nearer, F to move farther. With d, the second argument is an absolute position, such as is returned by 2POS?.

The third and fifth arguments set the acceleration at the beginning and deceleration at the end of the z

motion, though I have not yet figured out in what units. The fourth argument is the speed at which to move in units of a tenth of a micron per second (this is inconsistent with the distance measurements in the chassis!).

2POS? and 2MOV are inconsistent with the rest of the commands: the chassis responds to 2POS and 2MOV? as errors. If you issue 2MOV when the chassis is already executing another 2MOV command it returns 2MOV !,E02110. To stop a motion in progress, send 2STOP. It will return 2STOP +, and the currently executing motion will return an error with 2MOV !,E02133.

Autofocus

I only deal with the ZDC hardware autofocus system. Software autofocus is entirely independent of the microscope, and is inevitably a bad hack to work around the absence of a hardware autofocus. *All motorized microscopes should have a hardware autofocus.* Unfortunately, most microscope software persists in treating it like software autofocus.

ZDC works by firing a laser through the objective at the sample. The laser reflects off the near side of the coverslip (its interface with air) for air immersion and low magnification objectives, or off the far side of the coverslip where it meets the sample in oil and water immersion and high magnification objectives. One side effect of this is that the microscope knows which way it needs to move to go towards focus. This is distinctly different from software autofocus.

ZDC commands begin with 2AF. To set the near and far limits of the range in which the ZDC can search, use

```
2AFFLMT 536531
2AFNLMT 541531
```

You can give yourself plenty of space. Unlike software autofocus, hardware autofocus is slowed only by how far it is from correct focus, not by how large a search range you allow it. You also have to set the 2NEARLMT and 2FARLMT commands or the ZDC will fail with 2AF !,E02311.

At least before the first time using the ZDC in a session (though you can do so as often as you want), issue the magic command

```
2aftim 4
```

I couldn't find any conditions on my microscope which changed this, which makes it impossible to reverse engineer. Then you have to set up your objective with

```
2AFTBL nn
```

where nn is a number corresponding to your objective (see the table in the appendix for the command 2AFTBL for its values). Finally, to actually move the objective where it is focused on the edge of the coverslip, use

```
2AF SHOT
```

If it succeeds, it will reply 2AF +. If it fails, it returns an error code:

- 2AF !,E02312 = the coverslip's boundary is outside the far range limit

- 2AF !,E02313 = the coverslip's boundary is outside the near range limit
- 2AF !,E02331 = the ZDC could not find a coverslip boundary

The first two mean you are probably using too small a search range. If you get the third, or you get one of the first two and you know your search range is adequate, then make sure you have adequate oil with no bubbles on your objective, and that your sample is sufficiently rigid and heavy not to be shifted by the microscope. For example, ZDC usually fails if you are using only coverslips for your sample. If you mount your sample between a coverslip and a heavy glass slide instead of two coverslips, this often fixed the problem. Similarly, a colleague found that if she made her microfluidic devices much thicker (and so much heavier), the ZDC magically started working consistently.

Unknown

1SND0B 1SW 2AFSTS

Appendix: Command Reference

1BPORT

Switches to the back port of the chassis, if it has one. Not tested.

1CD

Usage: 1CD *pos*

Switch to position *pos* of the motorized condensor.

1LED

Usage: 1LED (1|0)

Open or close the external shutter for the fluorescent light source. Whether 1 corresponds to open or closed depends on your shutter. 1 corresponds to high voltage on the control line.

1LMP

Usage: 1LMP *inten*

Set the intensity of the transmitted light source to *inten*, measured in tenths of a voltage between 0V (1LMP 0) and 12V (1LMP 120).

1LMPSEL

Usage: 1LMPSEL *lamp*

Select illumination source *lamp*. Almost always *lamp* will be DIA.

1LMPSW

Usage: 1LMPSW (ON|OFF) Turn the currently selected lamp (usually the transmitted light source) on and off.

1LOG

Usage: 1LOG (IN|OUT)

IN means place all commands beginning with 1 under computer control, locking out the microscope chassis's physical controls unless specifically reenabled. OUT means to return control to the chassis's physical controls, and prevent the computer from making changes to the microscope's state. Mnemonic: "log in"

1MU

Usage: 1MU *pos*

Sets the filter wheel to position *pos*.

1OB

Usage: 1OB *pos*

Sets the objective nosepiece to position *pos*.

1PRISM

Usage: 1PRISM (1|2)

1 diverts the light path to the ocular. 2 diverts it to the camera. You should probably close the external fluorescence shutter with 1LED when you have diverted light to the ocular.

1SHUT1

Usage: 1SHUT1 (IN|OUT)

Opens (OUT) and closes (IN) shutter 1 (connected to port SHA1 on the chassis -- usually transmitted light).

1SHUT2

Usage: 1SHUT2 (IN|OUT) Opens (OUT) and closes (IN) shutter 2 (connected to port SHA2 on the chassis -- usually transmitted light).

1SNOB

Unknown.

1SW

Usage: 1SW (ON|OFF) Turns on/off using the buttons on the microscope chassis. When on, all button releases result in 1SW 0 being sent to the computer. Button presses result in 1SW *n*, where *n* is an integer indicating the button.

1UNIT

Usage: 1UNIT?

Returns a string describing the microscope. Our scope returns 1UNIT IX2,FRM,RV1,FO,MU6,HS.

1peekb

Usage: 1peekb *reg*

Looks at the value in register *reg*. CellR sends 1peekb D0003 during startup and receives the response 1peekb C7.

2AF

Usage: 2AF SHOT

Autofocus with the ZDC. For air immersion, low magnification objectives, it leaves the objective's focal plane at the air/coverslip interface nearest the objective. For water or oil immersion objectives, it leaves the focal plane at the interface farthest from the objective. If it succeeds in doing so, it returns 2AF +. If it fails, it returns one of the three following errors:

- 2AF !,E02312 = the coverslip's boundary is outside the far range limit
- 2AF !,E02313 = the coverslip's boundary is outside the near range limit
- 2AF !,E02331 = the ZDC could not find a coverslip boundary
- 2AF !,E02311 = you haven't set 2NEARLMT and 2FARLMT

Errors 2312 (outside the far range limit) and 2313 (outside the near range limit) can be fixed by extending the range with 2AFFLMT and 2AFNLMT. 2331 (could not find boundary) usually results from one of three things: your oil on the objective is inadequate or full of bubbles; your sample beyond the coverslip is dried out; or your assembly around your sample isn't heavy enough, such as if you put your sample between two 24x60mm coverslips instead of a coverslip and a heavy glass slide.

Before you execute 2AF SHOT, you need to send the commands 2AFFLMT, 2AFNLMT, 2AFTBL, and 2aftim.

2AFFLMT

Usage: 2AFFLMT *pos*

Set the far limit of the range of z-positions in which the ZDC can search for the coverslip position to *pos*, measured in hundredths of a micron, and from the same origin as the positions returns by 2POS?.

2AFNLMT

Usage: 2AFNLMT *pos*

Set the near limit of the range of z-positions in which the ZDC can search for the coverslip position to *pos*, measured in hundredths of a micron, and from the same origin as the positions returns by 2POS?.

2AFSTS

Unknown.

2AFTBL

Usage: 2AFTBL *num*

This configures the ZDC to use a particular objective. It must be done before using 2AF SHOT. The numbers for the different objectives are

<i>num</i>	Objective Type
75	LCPlanFl 20x
76	LCPlanFl 40x
79	LCPlanFl 60x
78	LUCPlanFl 40x
75	LUCPlanFIN 20x
76	LUCPlanFIN 40x
79	LUCPlanFIN 60x
31	PlanApo 100x O3
70	PlanApo 60x O TIRFM-SP
62	PlanApo 60x O/LSM
30	PlanApo 60x O3
60	PlanApo 40x W/LSM
30	PlanApoN 60x O
77	SLCPlanFl 40x
45	UApo 30x 3/340
47	UApo 40x 3/340
48	UApo 40x W
57	UPlanFl 60x OI3
36	UPlanAPO 20x
43	UPlanApo 100x OI3
39	UPlanApo 40x OI3
40	UPlanApo 60x
41	UPlanApo 60x W
42	UPlanApo 60x W/IR
58	UPlanFl 100x O3
55	UPlanFl 20x
71	UPlanFl 40x O-SP
56	UPlanFl 40x
58	UPlanFIN 100x O
58	UPlanFIN 100x OI
55	UPlanFIN 20x
39	UPlanFIN 40x O
40	UPlanFIN 60x
57	UPlanFIN 60x OI
43	UPlanSApo 100X O

36	UPlanSApo 20x
47	UPlanSApo 40x
30	UPlanSApo 60x O
41	UPlanSApo 60x W
56	UPlanFIN 40x

2aftim

Usage: 2aftim *n*

This is a magic command which I don't understand, but it needs to be sent to the ZDC before using 2AF SHOT. *n* is 4 as sent by CellM, and can be any of 1, 2, 3, or 4. Perhaps it's a timeout command for when the ZDC fails?

2FARLMT

Usage: 2FARLMT *pos*

Sets the farthest z position to which you can move the objective to *pos*, measured in hundredths of a micron from the same origin as 2POS?.

2JOG

Usage: 2JOG (ON|OFF) Enables (ON) and disables (OFF) the focus wheel on the microscope chassis after you have executed 2LOG IN. See 2JOGSNS to set the rate at which the wheel moves the objective.

2JOGSNS

Usage: 2JOGSNS *sens*

2JOGSNS sets the rate at which the focus wheel on the microscope chassis moves the objective nosepiece in the z direction. *sens* is an integer, at least 0 and no more than 10. The speed increasing exponentially with *sens*, from about 0.2um per turn of the wheel at *sens* = 0 to 800um per turn of the wheel at *sens* = 10.

2joglmt

Usage: 2joglmt (ON|OFF) If the chassis focus wheel is enabled (with 2JOG), ON prevents you from going beyond the limits set with 2NEARLMT and 2FARLMT with the focus wheel. OFF removes that limitation.

2LOG

Usage: 2LOG (IN|OUT)

IN means place all commands beginning with 2 under computer control, locking out the microscope chassis's physical controls unless specifically reenabled. OUT means to return control to the chassis's physical controls, and prevent the computer from making changes to the microscope's state. Mnemonic: "log in"

2MOV

Usage: `2MOV (N,F,d),x,start,speed,end` Moves the objective's z-position nearer (*N* for the first argument), farther (*F* for the first argument), or to an absolute position (*d* for the first argument). *x* is the distance to move or the absolute position to move to. *speed* is the speed at which to move in tenths (not hundredths!) of a micron per second. *start* and *end* control how much time at the beginning and end of the move is spend accelerating. Typically *start* is small, say 1, bringing the objective immediately up to speed, and *end* is larger, around 50, in hopes of getting better accuracy in the z motion. The four important error conditions for `2MOV` are:

- `2MOV !,E02414` - the motion reached the near limit (as set by `2NEARLMT`)
- `2MOV !,E02412` - the motion reached the far limit (as set by `2FARLMT`)
- `2MOV !,E02120` - the arguments to `2MOV` are invalid (no motion takes place)
- `2MOV !,E02133` - `2MOV` was ended by `2STOP`.

2NEARLMT

Usage: `2NEARLMT pos`

Sets the nearest z position to which you can move the objective to *pos*, measured in hundredths of a micron from the same origin as `2POS?`.

2POS?

Usage: `2POS?`

The microscope responds, `2POS pos`, where *pos* is an absolute z position of the objective measured in hundredths of a micron.

2STOP

Usage: `2STOP`

Stop any currently executing motion. Always responds with `2STOP +`. If there is a `2MOV` command in progress, it also aborts and returns an error condition with `2MOV !,E02133`.